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☰ < 100 % [Print] [Info] [List]

Safety Setting

Performance Safety Settings | **Human - Machine Safety Settings** | Safety IO Settings

Body Region Risk Setting X More Limit Setting X

1. Please set body regions that could be contacted by the robot in the collaborative workspace

Face & Head
 Neck Muscle
 Shoulder Joint and Muscle
 Upper Arm & Elbow Joint
 Lower Arm & Wrist Joint
 Lower Leg
 Thigh & Knee
 Hand & Fingers
 Abdominal Muscle
 Chest

2. Result

When robot enters the collaborative workspace, the path motion set with 100% speed will be automatically changed into mm/sec

When robot enters the collaborative workspace, the PTP motion set with 100% speed will be automatically changed into %

When robot enters the collaborative workspace, the speed setting will be achieved within ms

Enable G-Sensor ?

Input boxes highlighted with blue borders have been modified to be higher than the default values. Minimum possible contact area calculation have been removed. In this case, it's user's responsibility and liability to perform risk assessment and validate safety parameters due to the modification

This feature is designed for user to quickly set up an initial robot application in collaborative workspace following the biomechanical limits of each body region listed in ISO/TS 15066*. User should still perform risk assessment on real application before deployment. User shall take responsibility for human body region which is not

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


Plug & Play Options

- Grippers
 - Robotiq (commercialized by Omron)
 - OnRobot
 - Schmalz
 - Schunk
- Connectivité
 - ADLINK
 - Advantech
- Applications
 - PickIt: 3D Bin Picking
 - Kilews: Screw Driving
 - Basler: Vision
 - Force Sensor: Robotiq and ATI
- Autres
 - Murrplastik

Mobile Working Station (Uniquement TM5) **OMRON**

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

Poids: 133kg

IHM Tablette

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- Main Features:
 - 10.1" 16:10 WXGA 1280x800 500 nits LED Panel
 - 5 points projected capacitive touch
 - Full IP65 protection
 - Two external USB2.0 for data backup
 - Anti-vibration/shock IEC/EN 61131-2 compliance
 - EMC (IEC/EN 61000-6-2/4, IEC/EN 61131-2) compliance
 - Weight: ~1.5 kg without external cable

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Normes et Standards:

ISO 10218-1:2011

Robots and Robotic Devices - Safety requirements for industrial robots
(quasi machine)

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ISO 10218-2: 2011

Robots and robotic devices – Safety requirements for industrial robots. Robot systems and integration (machine)

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ISO / TS 15066

Technical Specification. Robots and Robotic devices – Collaborative Robots

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ISO 13849-1:2015

defines how a safety-related control system must be designed to achieve a specific performance level (Calcul des niveaux de performance)

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Safety Function Number	Name	Robot Stopping Function	Structure Category	PL
SF0	Robot Stick ESTOP	Emergency Stop	Cat. 3	d
SF1	User Connected ESTOP Input	Emergency Stop	Cat. 3	d
SF2	Encoder Standstill	Protective Stop	Cat. 2	d
SF3	User Connected External Safeguard Input	Protective Stop	Cat. 2	d
SF4	Joint Torque Monitoring	Protective Stop	Cat. 2	d
SF5	Joint Position Limit	Protective Stop	Cat. 2	d
SF6	Joint Speed Limit	Protective Stop	Cat. 2	d
SF7	TCP Speed Limit	Protective Stop	Cat. 2	d
SF8	TCP Force Limit	Protective Stop	Cat. 2	d
SF9	User Connected External Safeguard Input for Human-Machine Safety Settings	--	Cat. 2	d
SF10	Robot ESTOP Output	--	Cat. 3	d
SF11	User Connected External Safeguard Output	--	Cat. 3	d
SF12	User Connected External Safeguard Output for Human-Machine Safety Settings	--	Cat. 3	d
SF13	Robot Internal Protective Stop Output	--	Cat. 3	d
SF14	Robot Encoder Standstill Output	--	Cat. 3	d
SF15	User Connected Enabling Device Input	Protective Stop	Cat. 2	d
SF16	User Connected ESTOP Input without Robot ESTOP Output	Emergency Stop	Cat. 3	d
SF17	Cartesian Limit	Protective Stop	Cat. 2	d

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TM Robot

ISO/TS 15066 oriented design

1. Check risk parts
2. Robot will adjust and set the speed limit automatically

Other Robots

Each of the following joint limits can be configured independently.					
JOINTS	Range	Normal Mode		Reduced Mode	
		Minimum	Maximum	Minimum	Maximum
Base	363° - 363°	363°	363°	363°	363°
Shoulder	363° - 363°	363°	363°	363°	363°
Elbow	363° - 363°	363°	363°	363°	363°
Wrist 1	363° - 363°	363°	363°	363°	363°
Wrist 2	363° - 363°	363°	363°	363°	363°
Wrist 3	363° - 363°	363°	363°	363°	363°

Body region	Specific body area	Quasi-static contact		Transient contact	
		Maximum permissible pressure* Fv N/cm²	Maximum permissible force* Fv N	Maximum permissible pressure multiplier† Fv	Maximum permissible force multiplier† Fv
Skull and forehead	1. Middle of forehead	130	130	not applicable	not applicable
	2. Temple	110	110	not applicable	not applicable
Face*	3. Mastoid muscle	110	65	not applicable	not applicable
	4. Neck muscle	140	150	2	2
Neck	5. Seventh neck muscle	210	150	2	2
	6. Shoulder joint	160	210	2	2
Back and chest	7. Fifth lumbar vertebra	210	210	2	2
	8. Sternum	120	140	2	2
Arm	9. Pectoral muscle	170	110	2	2
	10. Abdominal muscle	140	110	2	2
Pelvis	11. Pelvic bone	210	180	2	2
	12. Deltoid muscle	190	150	2	2
Anatomical joint	13. Humerus	220	150	2	2

1. Understand all parameter, their physical meanings, and the expected response mechanism.
2. Read and understand the detail knowledge of related international standards, make clear the regulations and domain knowledge of human-robot collaboration safety.
3. Become an expert of the above two domain.
4. Trial and error for the safety parameters.
5. Design and arrange experiments, e.g. using force sensors, load cells...to test if shock/clamping force/pressure are matching with safety regulations.
6. Or just try the pain by your self.
7. Congratulation! You finally set the robot and can publish a paper for it!

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Safety IO Settings

Cartesian Limit

Enable Cartesian Limit and Disable Operation Space ?

Limits of position

TCP elbow

Cartesian Setting

Cube

x upper bound 500 x lower bound -500

y upper bound 3000 y lower bound -3000

z upper bound 3000 z lower bound -3000

θ_z 0

Cylinder


Default

Save

Applications

Domaines d'application

OMRON

Pick & Place

Assemblage

Machine Tending

Inspection

Palettisation

Vissage

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MoMa : Mobile Manipulator